

LeakControl

Automated Endoscope Air Leak Tester

The only **air leak tester** that finds leaks early on when **repairs are minor and less costly**, minimizing the level of damage caused by submersion.

Prevent up to 90% of flexible endoscope costs



Comply to
regulations
and IFU



Eliminate the
risk of cross
infection



Prevent 90% of
flexible endoscope
costs



No more large repairs
on your flexible
endoscopes



No more
surprises



**What can Automated Endoscope
Leak Testing do for you?**



What will Automated Air Leak Testing do for you?

No more large repairs

With LeakControl as part of your manual pre-cleaning procedure, it prevents water ingress in your scopes and the large and costly repairs that result from it.

Eliminating cross infection

Because LeakControl detects the leaks that would otherwise go unnoticed, it prevents organic materials from being able to enter the scope or build up around the leak.

This results in the prevention of cross infection between patients.

Crucial insights

All test data is stored in our cloud platform EndoscopeManager. This ensures regulatory compliance and gives you crucial insights in your inventory, allowing you to track your light guides performance over time.


Overall benefits

- **Eliminate the risk of cross infection**
- **Reduce flexible endoscope costs by 90%**
- **Automatically comply with regulations**
- **Uninterrupted disinfection process**
- **Increase patient safety**
- **No more surprises**

What does LeakControl Test?

LeakControl tests your endoscopes by measuring pressure over time. On your device it will simply tell you if your scope has a leak.

In EndoscopeManager you can see the pressure over time, and the pressure drop in percentage.



LeakControl connected to EndoscopeManager.

Quality control, Predictability and risk management.

- **Accurate traceability** of scopes and parameters:
Historical data
- Increases **Predictability, Traceability and Quality Control**
- Optional **automated workflows**:
Notifications and/or autostarts and follow ups (e.g. repair/replace)
- **Improves logistics**, repair, replacement & delivery process and can be used by Hospitals, Manufacturers & Repair shops
- **Stores & evaluates results** from LeakControl

