




ScopeControl[®]

100% quality guaranteed in 3 minutes

The only fully automated rigid endoscope testing system that ensures your scopes are in perfect shape and ready for OR use.



What does ScopeControl[®] measure?

 *Field of View*

 *View Angle*


 *Light Fibres*

 *Light Transmission*

 *Colour Correctness*

 *Focus*

• *Easy to learn!* • *Fast to use!* • *100% quality assured!*



EndoscopeManager benefits:

- *Insights in production*
- *Clear report function*
- *Check repairs*
- *Follow lifetime and costs*



What will ScopeControl® do for me?

- ✓ Reduce OR costs through improved efficiency
- ✓ Compare cost of device usage per operation
- ✓ Identify minor issues before they become major repairs
- ✓ Reduced delays in OR and happy surgeons
- ✓ Checks your scopes daily before use
- ✓ Better patient safety and reduced infection risk
- ✓ Track inventory and scope performance over time
- ✓ Reduce time testing scopes and eliminate guesswork



ScopeControl®

Testing in Three Easy Steps...

- 1) Place the endoscope inside ScopeControl and close the lid
- 2) Use the touch screen or a scanner to identify the scope
- 3) Press start




View a demo of the system, scan code:



In 1-3 minutes, the system will produce a simple pass fail result from the test

Test results and lens pictures are automatically saved

Management reports and remote access available from any place, any time



For further information or to find your distributor please visit our website or contact us:



DOVIDEQ medical BV | Keulenstraat 8G 7418 ET Deventer NL | KVK 52693031 | VAT NL8505.55.279.B01
PHONE +31 (0) 570760800 | FAX +31 (0) 570760809 | E-MAIL info@dovideqmedical.com | WEB www.dovideqmedical.com

No liability is accepted for errors and misprints. All specific product conditions stated in this document are subject to minor technical modifications and design deviations.